



Product Specifications

Item No.	Description	A	B	C	D	E
		Overall Length (mm)	Diameter of tip in middle phalanx (mm)	Diameter of tip in proximal phalanx (mm)	Length of angled implant in middle phalanx (mm)	Length of angled implant in proximal phalanx (mm)
61-253519	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 19mm Straight	19	3.5	2.5	-	-
61-253522	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 22mm Straight	22	3.5	2.5	-	-
61-253525	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 25mm Straight	25	3.5	2.5	-	-
61-253530	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 30mm Straight	30	3.5	2.5	-	-
61-253535	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 35mm Straight	35	3.5	2.5	-	-
67-253519	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 19mm Angled	19	3.5	2.5	8	11
67-253522	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 22mm Angled	22	3.5	2.5	9	13
67-253525	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 25mm Angled	25	3.5	2.5	9	16
67-253530	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 30mm Angled	30	3.5	2.5	9	21
67-253535	ARROW-LOK Hybrid: Dia. 2.5-3.5mm / Length 35mm Angled	35	3.5	2.5	9	26
61-3513	ARROW-LOK: Dia. 3.5mm / Length 13mm Straight	13	3.5	3.5	-	-
61-3516	ARROW-LOK: Dia. 3.5mm / Length 16mm Straight	16	3.5	3.5	-	-
61-3519	ARROW-LOK: Dia. 3.5mm / Length 19mm Straight	19	3.5	3.5	-	-
61-3522	ARROW-LOK: Dia. 3.5mm / Length 22mm Straight	22	3.5	3.5	-	-
67-3516	ARROW-LOK: Dia. 3.5mm / Length 16mm Angled	16	3.5	3.5	6	10
67-3519	ARROW-LOK: Dia. 3.5mm / Length 19mm Angled	19	3.5	3.5	8	11
67-3522	ARROW-LOK: Dia. 3.5mm / Length 22mm Angled	22	3.5	3.5	9	13

*This illustration represents the ARROW-LOK design with 3.5mm diameter arrowhead tips at each end. The ARROW-LOK Hybrid design features a 3.5mm diameter arrowhead tip at the end labeled "D" which is inserted in the middle phalanx and a 2.5mm diameter arrowhead tip at the end labeled "E" inserted in the proximal phalanx.